

Barchester Hospital/Childrens Block

Period: 12 months up to and including 07/2011

Level 1

Room	Source type	Ref No	System No	Temp Cat	Task	Aug	Sept	Oct	Nov	Dec	Jan	Feb	Mar	Apr	May	June	July
A1/001 Boiler House F	Bib Tap	20394		C	Check cold temp	18.1,	18.9,	17.1,	13.5,		12.0,	11.0,	13.0,	13.0,	13.0,	15.9,	
						18.0,	20.0,	15.2,	11.9		12.0	11.0,	9.0	11.0,	14.0,	17.6,	
						18.9	19.2,	12.9			13.0		12.1,	14.0	19.1,		
							18.4,						11.0			18.0,	
							17.9									19.4	

Level 2

Room	Source type	Ref No	System No	Temp Cat	Task	Aug	Sept	Oct	Nov	Dec	Jan	Feb	Mar	Apr	May	June	July
A2/002 General Paed Lab	Cold Tap	20086		C	Check temp (cold)	20.0	20.0	19.0			13.6	9.0	13.0	12.0	16.0	18.0	
	Hot Tap	20079		H	Check temp (hot)	57.0	60.0	60.5			64.3	65.0	66.0	66.0	65.0	65.0	
A2/003 Lab Office	Cold Tap	20088		C	Check temp (cold)	19.1		19.0			18.0	9.0	22.0*				
	Hot Tap	20087		H	Check temp (hot)	47.1*		52.0			54.0	63.0	43.0*				
A2/019 Metabolic Lab	Cold Tap	20090		C	Check temp (cold)	19.1	20.0	20.0			13.0	11.0	17.0			17.0	
	Hot Tap	20089		H	Check temp (hot)	54.1	59.0	52.0			64.4	54.0	65.0			65.0	
A2/102 Domestics Store	Cold Tap	20092		C	Check temp (cold)	20.0*	21.0*	21.0*			19.0	14.0	10.0	18.0	20.0	24.0*	
	Hot Tap	20091		H	Check temp (hot)	59.1	53.0	58.0			58.0	56.0	59.0	63.0	64.0	65.0	

Level 3

Room	Source type	Ref No	System No	Temp Cat	Task	Aug	Sept	Oct	Nov	Dec	Jan	Feb	Mar	Apr	May	June	July
A3/001 Diagnostic Lab	Cold Tap	20094		C	Check temp (cold)	19.0	19.4	18.2			13.8	12.0	15.0	14.0	18.9	17.9	
	Hot Tap	20093		H	Check temp (hot)	56.0	63.7	62.4			65.1	65.0	66.0	65.0	65.0	65.4	
A3/005 Reg Haemostatis Lab 1	Cold Tap	20096		C	Check temp (cold)	19.0		18.7			16.0	10.9	13.0	12.0	20.0	17.0	
	Hot Tap	20095		H	Check temp (hot)	54.0		62.1			65.0	65.0	62.0	65.2	65.0	55.0	
A3/024 Molecular biology	Cold Tap	20098		C	Check temp (cold)	19.1	20.0	19.1			14.0	15.9	13.0	11.7	20.0	17.0	
	Hot Tap	20097		H	Check temp (hot)	54.0	63.0	60.3			63.5	55.0	64.5	65.4	64.0	65.0	

Level 4

Room	Source type	Ref No	System No	Temp Cat	Task	Aug	Sept	Oct	Nov	Dec	Jan	Feb	Mar	Apr	May	June	July
A4/007 Media Room	Cold Tap	20102		C	Check temp (cold)	22.0*	19.4	0.0			12.0	12.0	13.0	13.0	18.0	18.9	
	Hot Tap	20101		H	Check temp (hot)	56.0	65.0	60.7			65.0	64.0	65.0	65.2	66.0	65.4	